

## 1. **Applicable scope**

Can be used to treat intervertebral disc herniation.

## 2. **Contraindication**

Spinal spondylolisthesis and spinal instability;

Psychiatric patients;

Coagulation disorders, severe heart and lung diseases, and pregnancy.

## 3. **Warnings and Cautions**

1. Pay attention to check the device packaging before use, if the packaging is damaged or sterilization expired, do not use.

2. This tool is disposable, and disposable tools can cause environmental pollution. It should be disposed of in accordance with relevant national or industry regulations to avoid environmental pollution.

3. This product is only for one-time use. Repeated sterilization poses a risk of equipment damage or cross infection.

4. Only the trained/experienced physicians should use this device.

5. The device can not cause death of any user or serious deterioration in health status.

6. This instruction manual must be read carefully before use.

7. When disposing of medical waste, proper protection should be taken to avoid infection and microbial hazards, and attention should be paid to avoiding injury from sharp objects such as Puncture Needle, Trephine.

## 4. **Storage**

The sterilized Percutaneous Discectomy Instruments Pack should be stored in a well-ventilated room with 0-40°C , relative humidity not more than 80%, free of corrosive gas.

## 5. **Complications/Side Effects**

1. Postoperative Dysesthesia
2. nerve damage or paralysis, weakness
3. infection

4. worse pain

5. bleeding

## 6. Labels and Symbols:



Medical device



Model number



Manufacturer



Date of Manufacture




Use-by date



Batch Code



Sterilized using ethylene oxide



Double sterile barrier system



Do not re-use



Do not use if package is damaged



Temperature limit




Caution



Humidity limitation



Do not resterilize



**Company Name : Lotus NL B.V.**

**Company Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherl**

Date: Sep.,14,2023