

AESTHETIC & MEDICAL

Leading IPL, Laser, Slimming Machine Manufacturer

6D Lipo Laser



sunrise



6D lipo laser body slimming beauty machine




How does it work

6D Lipo Laser

635/650nm laser wavelength has been proven to shrink fat cells by tricking them into making them think that you are working out. This wave will penetrate 8-10mm to the fat layers, the cells release stored fat in order for your body to burn it as energy.

780/808nm laser wavelength penetrate 10-13mm deeper to the fat layers, making the ablation fat and toxins eliminating from the body quickly, and strengthen the body shaping effect.

940/980nm laser wavelength penetrate 13mm-16mm deeper to the fat layers, burn out the stubborn cellulite into triglyceride then dissolve it, accelerate lymphatic drainage,



Advantages


- 1.** Lipo laser machines are 100% non-invasive, laser based, spot fat reduction and body contouring systems.
- 2.** New design, its wavelengths 650nm & 780nm & 808nm & 940nm & 980nm can work together or independently.
- 3.** 635nm X 5 lights + 650nm X 5 lights + 780nm X 5 lights + 808nm X 4 lights + 940nm X 5 lights + 980nm X 4 lights, every big pad with 28 lamps
- 4.** All the diode lights are imported from Japan Mitsubishi, the power is 209MW, it is highest in china market.
- 5.** Germany imported high quality electronic board, every pad has its own power source and air cooling system, which is convenient to maintain, and prolongs lifespan of diode lasers.
- 6.** The pads can be connected tightly by bandage, which is more effective to lose weight.

Specification


6D Lipo Laser			
Wavelength	635nm、 650nm、 780nm、 810nm、 940nm、 980nm		
Diode Power	Single diode 100-200mw		
Light No.	28	Screen	12.1Inch
Pads	4,6,8,10,12,14		
Machine size	55*30*43	Packing size	66*51*51cm
Input	AC110V-230V	Power	500VA
Working mode	continours / pulse setting		




**210mv for each diode laser light.
A pad with 28 diode laser lights will be
5880mv in total.**



Screen Show




Pads show




635nmX5 lights
650nmX5 lights
780nmX5 lights
808nmX4 lights 940nmX5 lights
980nmX4 lights every big pad
with 28 lamps

Treatment Procedure

- ▲ Before any treatments please refrain from eating a large meal and avoid caffeine and alcohol at least 3 hours before.
- ▲ Shower before treatments and make sure no cream or perfumed products have been applied to the skin for at least 24 hours before.
- ▲ Don't need to put any of gel or essential oil onto body , just need to be sure that parameter is adjust to the suitable level for different client.
- ▲ Suggest using 940nm first for fat dissolving , then 650nm for consolidating.



Treatment




Before and After



Before

4 Weeks After

6 Weeks After



Before

After

Before

After



Before

After

Before

After

THANKS

>Web: www.beautysunrise.com

>ADD 1: JinRi Science & Technology zone,Daxing District,Beijing,China

>ADD 2: Room 803, unit 1, buliding 40,Zaoyuan Dongli, Daxing District Beijing,China

Contcat Us

