

PEEK Interference Screw System 聚醚醚酮钉鞘固定系统

**Excellent Biocompatibility
Fixation Strength**

Comprehensive Healing Support

High

聚醚醚酮（PEEK）材料卓越的生物相容性，不显影，强度高且弹性模量与骨骼接近，安全可靠

PEEK (polyether ether ketone) offers excellent biocompatibility, radiolucency, high mechanical strength, and an elastic modulus similar to bone, ensuring safety and reliability.

系统采用模块化设计，界面螺钉既可单独使用，也可与外鞘配合使用，为外科医生提供了灵活的术中选择

The system features a modular design. The interference screw can be used independently or in combination with an outer sheath, providing surgeons with flexible intraoperative options.

外鞘尾部导向螺纹设计，螺钉植入方便，固定更牢固

The sheath features a threaded trailing end guiding design, facilitating easier screw insertion and improving fixation stability.

界面螺钉与外鞘配合，避免螺钉对移植物的切割损伤，360°最大化隧道内移植物接触，提供更为可靠的固定效果

The interference screw works in conjunction with the outer sheath to prevent graft laceration caused by screw insertion. It maximizes 360° graft-to-tunnel contact, providing more reliable fixation.