



微界
V Plus




Electronic Endoscopic Imaging System

Main Machine



- The main machine is equipped with a button board to adjust the working state of the machine.
- Support high-definition output, clear observation of the surgical process.
- The Main machine has the function of saving, recording and exporting the images.

Endoscopic end Angle classification



0°endoscope



30°endoscope

Disposable Electronic arthroscope



- The field of view is wide and clearer.
- The part touched human body will not be over-temperature, which is safe and reliable.
- Ergonomically designed handle to adapt to various holding habits.

Model

Specification(mm)

JWAA001	Φ4.0	0°/100°/140°	60-280 (interval 5)
JWAA002	Φ5.2	30°/100°/140°	60-280 (interval 5)




关注微界公众号



V PLUS business email: jxwj@weijie.co

Disposable Electronic spinal endoscope



- The field of view is wide and clearer.
- The part touched human body will not be over-temperature, which is safe and reliable.
- Ergonomic design handle, more comfortable to hold, more easy to use.

Model	Specification(mm)		
JWAW001	$\Phi 6.3 \times \Phi 3.7$	23°/120°	180
JWAW002	$\Phi 5.8 \times \Phi 3.1$	15°/120°	125

Open channel A



- In minimally invasive spine surgery, the supporting camera system cooperates with the orthopedic retractor, and the appropriate size of the opening channel A is selected to be inserted into the incorporation, and the surgical working channel is established.
- Open channel A is available in 4 sizes to suit different patient. The channel scope can be inserted at the same time to cooperate with the operation, with sufficient illumination and a clearer surgical vision.

Model:FAU/Q004

Specification(Length):	40/50/60/70 mm
Diameter:	$\Phi 21$ mm

Disposable Electronic channel endoscope



- Use with open channel A.
- The field of view is wide and clearer.
- The part touched human body will not be over-temperature, which is safe and reliable.

Model	Specification(mm)		
JWAP005	$\Phi 4.0$	0°/100°	30

Arthroscopic Cannulated



Model:JWAT002 Specification(mm)

$\Phi 4.0 \times 110 \times 0^\circ$
$\Phi 4.0 \times 142 \times 0^\circ$
$\Phi 4.0 \times 150 \times 0^\circ$
$\Phi 4.0 \times 175 \times 0^\circ$
$\Phi 4.0 \times 280 \times 0^\circ$
$\Phi 5.2 \times 110 \times 30^\circ$
$\Phi 5.2 \times 142 \times 30^\circ$
$\Phi 5.2 \times 150 \times 30^\circ$
$\Phi 5.2 \times 175 \times 30^\circ$
$\Phi 5.2 \times 280 \times 30^\circ$
$\Phi 6.0 \times 150 \times 0^\circ$
$\Phi 6.0 \times 175 \times 0^\circ$

It is used in conjunction with the endoscope to isolate the endoscope from the patient during surgical approach or examination to prevent cross-infection.

The above product specifications are partial specifications.
If you have special requirements, please contact us.

V PLUS business email:jxwj@weijie.co



关注微界公众号