

Cryotherapy Chamber

(Electric Style)

Model:uCyro E1



About Cryotherapy Therapy :

Ultra-low temperature therapy technology is a secret weapon for anti-aging, which narrows blood vessels and slows down blood flow and metabolism through the extreme low temperature inside the cabin, achieving anti-aging effects.

Due to the body cooling caused by whole body cryotherapy. it leads to an increase in basal metabolic rate, single dose whole body cryotherapy has been proven to burn an additional 500-800 calories in the body, The base metabolic rate measured by indirect calorimetry shows that even within a few weeks after a course of treatment, it can increase by 200-300 calories per day. Without changing diet, this means that you can burn 2 pounds of body fat per month without dieting or exercising.

Features :

- Automatic shutdown at night
- Can set pre cooling for startup
- 70°C or Other Minimum Temperature Is Customizable
- with secondary interface
- Both interfaces are adjustable

Specification:

Index	Content
Model Name	uCyro E1
Chamber Style	Enclosed Style
Size	1450(W)*1000(D)*2400(H)
Weight Capacity	380KG
Color	Default black or white (customizable)
Material	stainless steel
Power(W)	≅ 8500W
System	One Click Start/Stop
Control Panel	with secondary interface Both interfaces are adjustable
Intelligent Control	1. Automatic shutdown at night 2. Can set pre cooling for startup
Minimum Temperature	-70℃ Or Other Minimum Temperature Is Customizable
Single Use Time	180s
Operation System	Intelligent Operating System One Click Start/Stop
Time Setting Range	180s
Temperature Display	From room temperature to -70℃
Screen Size	10inch
Power supply voltage	110V-220V/50Hz
Power(W)	≅ 8500W

Configuration List

Index	Specification	Unit	QTY
Chamber	1450(W)*1000(D)*2400(H)	set	1
vacuum hose		set	1
power cable		set	1
Operation and Maintenance Manual		set	1

Chamber:



UI:





DR.HUGO

Your Therapy Expert

POWER BY:



LANNX Bio&Med Co.,Ltd

Website: www.lannx.net/cryotherapy

E-Mail: sale@lannx.net

-Search our channel in YouTube to view operation video:

<https://www.youtube.com/@cyrotherapychamber>