


Cryotherapy Chamber

(Portable Style)

Model:uCyro P8



About Cryotherapy Therapy:

Ultra-low temperature therapy technology is a secret weapon for anti-aging, which narrows blood vessels and slows down blood flow and metabolism through the extreme low temperature inside the cabin, achieving anti-aging effects.

Due to the body cooling caused by whole body cryotherapy, it leads to an increase in basal metabolic rate, single dose whole body cryotherapy has been proven to burn an additional 500-800 calories in the body. The base metabolic rate measured by indirect calorimetry shows that even within a few weeks after a course of treatment, it can increase by 200-300 calories per day. Without changing diet, this means that you can burn 2 pounds of body fat per month without dieting or exercising.

Features:

- The soft chamber and the control system are designed to be separated but can be stored together. It is hand-pulled, portable, and vehicle-mounted. The whole device weighs only 53kg and does not require manual fixation after inflation. It occupies only 1.4 square meters.
- Therapy Console: The height and screen tilt angle conform to ergonomics, and the handle is streamlined and integrated, with excellent touch feel. It can be placed sideways and transported in the trunk, making it highly portable.
- Maximum storage: The left side can store the foldable soft cabin, liquid nitrogen outlet pipe, liquid nitrogen hose, handheld emergency stop switch, manual, etc., and the right side can be opened for maintenance of the electric control cabinet.
- Portable soft chamber: Low temperature resistant, wear resistant, waterproof, lightweight; grid



reinforcement, frame support, no need for manual fixation, stable structure after inflation, easy to store after deflation; dedicated hidden zipper goes straight to the top of the soft cabin to maintain beauty and facilitate users to enter and exit the soft cabin without bending over.

-Dual scene use: For short-distance transportation, the therapy console (stored in one) can be transported by taxi. For long-distance transportation, it is equipped with a locking trolley aviation box, and the whole machine can be completely stored in the box for flight consignment and logistics vehicle transportation. It is truly portable and has a device positioning function to monitor the location of the device in real time (optional) automatic monitoring system.

-Safe to use: Dual temperature PT100 probe sensor, sensitive temperature measurement, provides dual low temperature protection, and has sound and light alarm function; emergency stop switch, if the user feels the temperature is too low or feels unwell during use, the user can press the emergency stop switch to stop the machine to ensure personal safety

-Personalization: The size, color and size of the soft cabin can be customized; the Velcro can be customized with different logos according to customer needs.

Specification:

Index	Content
Model Name	uCyro P8
Size	Chamber: Φ 1200mmx1600mm Control system: 400x600x1000mm
Weight Capacity	53KG
Color	blue
Maximum capacity	1 person
Temperature Range	-180~-90℃
Liquid nitrogen consumption	2L/Time
Cold therapy duration	3-5 minutes(two modes can be freely switched, and can be stopped at any time during the process)
Single liquid nitrogen consumption	2-3L
Pre cooling time	No need for pre-cooling, you can start the treatment directly at room temperature. Or start pre-cooling, the pre-cooling time is less than 90 seconds
Refrigerant	liquid nitrogen
Power supply voltage	110V-220V/50Hz
Thermal insulation coefficient	0.01w/(m.k)

Configuration List

Index	Specification	Unit	QTY
Chamber	Chamber: Φ 1200mmx1600mm	set	1
Control system	400x600x1000mm	set	1
vacuum hose		set	1
power cable		set	1
Operation and Maintenance Manual		set	1



DR.HUGO

Your Therapy Expert

POWER BY:



LANNX Bio&Med Co.,Ltd

Website: www.lannx.net/cryotherapy

E-Mail: sale@lannx.net

-Search our channel in YouTube to view operation video:

<https://www.youtube.com/@cryotherapychamber>