

Instruction for Use



Device Description and Specification

Product

Gauze Sponge without X-ray

Photographs



Model

Gauze Sponge without X-ray is containing one model with different specifications, and the company specifications are described as follows.” MG “is the abbreviation of the company, “GS “is the abbreviation of the Gauze Sponge without X-ray,

Items	Model	L	W	Layer
1	MG-GS	5	5	4P
2	MG-GS	5	5	8P
3	MG-GS	5	5	12P
4	MG-GS	5	5	16P
5	MG-GS	5	5	32P
6	MG-GS	75	5	4P
7	MG-GS	75	5	8P
8	MG-GS	75	5	12P
9	MG-GS	75	5	16P
10	MG-GS	75	5	32P
11	MG-GS	10	5	4P
12	MG-GS	10	5	8P
13	MG-GS	10	5	12P

14	MG-GS	10	5	16P
15	MG-GS	10	5	32P
16	MG-GS	15	5	4P
17	MG-GS	15	5	8P
18	MG-GS	15	5	12P
19	MG-GS	15	5	16P
20	MG-GS	15	5	32P
21	MG-GS	20	5	4P
22	MG-GS	20	5	8P
23	MG-GS	20	5	12P
24	MG-GS	20	5	16P
25	MG-GS	20	5	32P
26	MG-GS	6	6	4P
27	MG-GS	6	6	8P
28	MG-GS	6	6	12P
29	MG-GS	6	6	16P
30	MG-GS	6	6	32P
31	MG-GS	75	6	4P
32	MG-GS	75	6	8P
33	MG-GS	75	6	12P
34	MG-GS	75	6	16P
35	MG-GS	75	6	32P
36	MG-GS	8	6	4P
37	MG-GS	8	6	8P
38	MG-GS	8	6	12P
39	MG-GS	8	6	16P
40	MG-GS	8	6	32P
41	MG-GS	9	7	4P

42	MG-GS	9	7	8P
43	MG-GS	9	7	12P
44	MG-GS	9	7	16P
45	MG-GS	9	7	32P
46	MG-GS	10	7	4P
47	MG-GS	10	7	8P
48	MG-GS	10	7	12P
49	MG-GS	10	7	16P
50	MG-GS	10	7	32P
51	MG-GS	75	75	4P
52	MG-GS	75	75	8P
53	MG-GS	75	75	12P
54	MG-GS	75	75	16P
55	MG-GS	75	75	32P
56	MG-GS	10	75	4P
57	MG-GS	10	75	8P
58	MG-GS	10	75	12P
59	MG-GS	10	75	16P
60	MG-GS	10	75	32P
61	MG-GS	15	75	4P
62	MG-GS	15	75	8P
63	MG-GS	15	75	12P
64	MG-GS	15	75	16P
65	MG-GS	15	75	32P
66	MG-GS	20	75	4P
67	MG-GS	20	75	8P
68	MG-GS	20	75	12P
69	MG-GS	20	75	16P

70	MG-GS	20	75	32P
71	MG-GS	8	8	4P
72	MG-GS	8	8	8P
73	MG-GS	8	8	12P
74	MG-GS	8	8	16P
75	MG-GS	8	8	32P
76	MG-GS	10	8	4P
77	MG-GS	10	8	8P
78	MG-GS	10	8	12P
79	MG-GS	10	8	16P
80	MG-GS	0	83	2P
81	MG-GS	9	9	4P
82	MG-GS	9	9	8P
83	MG-GS	9	9	12P
84	MG-GS	9	9	16P
85	MG-GS	9	9	32P
86	MG-GS	10	9	4P
87	MG-GS	10	9	8P
88	MG-GS	10	9	12P
89	MG-GS	10	9	16P
90	MG-GS	10	9	32P
91	MG-GS	15	9	4P
92	MG-GS	15	9	8P
93	MG-GS	15	9	12P
94	MG-GS	15	9	16P
95	MG-GS	15	9	32P
96	MG-GS	10	10	4P
97	MG-GS	10	10	8P

98	MG-GS	10	10	12P
99	MG-GS	10	10	16P
100	MG-GS	10	10	32P
101	MG-GS	12	51	04P
102	MG-GS	12	51	08P
103	MG-GS	12	51	012P
104	MG-GS	12	51	016P
105	MG-GS	12	51	032P
106	MG-GS	13	10	4P
107	MG-GS	13	10	8P
108	MG-GS	13	10	12P
109	MG-GS	13	10	16P
110	MG-GS	13	10	32P
111	MG-GS	15	10	4P
112	MG-GS	15	10	8P
113	MG-GS	15	10	12P
114	MG-GS	15	10	16P
115	MG-GS	15	10	32P
116	MG-GS	20	10	4P
117	MG-GS	20	10	8P
118	MG-GS	20	10	12P
119	MG-GS	20	10	16P
120	MG-GS	20	10	32P

Intended Use

Gauze Sponge without X-ray is used to cover/clean wounds or abrasions and absorb their exudates, absorb body-surface exudates. It is only intended to be used as a mechanical barrier; and not contact mucous membrane or manage the micro-environment of injured skin.



Document No.: MG/CE-MDR-03-02-02 A/0

Intended Users

Medical professionals in a specialty, clinical nurses

Intended Patient Populations

Patients with exudate from wound, which including all age groups.

Age: no restrictions

Weight: no restrictions

Gender: no restrictions

Nationality: multiple

Intended using environment

Hospitals, health centers, clinics.

Side effect

Wound infection caused by improper operation.

Medical Conditions

With wound exudate

Indications

Patients with exudate from wound.

Contraindications

None

Not applicable

People allergic to cotton

Sterilization

EO sterilization

Shelf life:5 years

Cautions

- 1) This product is sterilized for single use only, discard after single use. Please do not repeat sterilization, do not use it once again.
- 2) Should be used the product immediately after opening the individual package, please destroy when it was used.
- 3) Do not use it if the inner package has been previously opened or damaged, or

package be affected with damp, or the product went moldy.

4) Do not use after expiry date. If over the expiry date, please discard it.

Warning

It is forbidden for people who are allergic to cotton yarn

Application method

- 1) Peel off the package. Place Gauze Sponge without X-ray on the wound to absorb the exudate from the wound
- 2) Discard in the trash after use.

Properties of Product:

It has a combination with non-toxic, medical, environmental protection features. It has the application of special technology, the unique process, through the multi-process processing so that the product mesh and elasticity shall always keep uniform, better permeability.

Transportation:

It should be protected from rain and moisture during transportation.

Storage

Store in a cool, dry and dark place. Avoid direct sunlight

Symbols meaning

Signs	Description	Signs	Description	Signs	Description
	Trademark of Manufacturer		Batch number		Date of Manufacture
	Medical device		Model number		Country of manufacturer
	Use-by date		Sterilized using ethylene oxide		Caution
	Non-pyrogenic		Catalogue number		Unique Device

Signs	Description	Signs	Description	Signs	Description
					Identifier
	Single sterile barrier system		Double sterile barrier system		Keep away from sunlight
	Do not reuse		Do not use if the package is damaged		Keep dry
	Do not resterilize		Authorized representative in the European Community		Manufacturer
	CE marking of conformity, and Notified Body Code		Consult instructions for use	/	/

Manufacturer

Name: Haian Medigauze Co., Ltd.

Address: Jiaoxie Industry Park, 226633 Haian, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-513-88240958

Web: www.medigauze.com

E-mail: produce@medigauze.com

SRN: CN-MF-000009960

European Authorized Representative

Name: Shanghai International Holding Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537 Hamburg, Germany

E-mail: shholding@hotmail.com



Document No.: MG/CE-MDR-03-02-02 A/0

SRN: DE-AR-000000001