



Instruction for Use

Product

Non-woven gauze pads

Specification

MG-NS-05054P, MG-NS-05056P, MG-NS-05058P, MG-NS-75754P, MG-NS-75756P,
MG-NS-75758P, MG-NS-10104P, MG-NS-10106P, MG-NS-10108P, MG-NS-10204P,
MG-NS-10206P, MG-NS-10208P

Basic UDI-DI

697313440MGNS008N6

Device figure



Materials

Viscose polyester fiber

Caution

- 1) This product is for single use only, please do not reuse it.
- 2) Please check the integrity of the package before use, if the package is damaged, please do not use it.
- 3) Please read the instructions carefully before use. If the product has expired, please do not use it.
- 4) Please inform the manufacturer and competent authority in case of any adverse events related to the equipment failure.

Storage conditions

The product shall be stored in a dry, ventilated, non-corrosive environment and kept away from fire and flammable materials.

Shelf-life

5 years

Sterility

EO sterilization

Intended User

Medical professionals in a specialty, clinical nurses

Intended Patient Populations

Patients with exudate from wound, which including all age groups.

Age: no restrictions

Weight: no restrictions

Gender: no restrictions

Nationality: multiple

Intended using environment

Hospitals, health centers, clinics.

Side effect

Wound infection caused by improper operation.

Medical Conditions

With wound exudate

Indications

Patients with exudate from wound.

Contraindications

None

Definitions of Signs

Signs	Description	Signs	Description	Signs	Description
	Trademark of Manufacturer		Batch number		Date of Manufacture
	Medical device		Model number		Country of manufacturer
	Use-by date		Sterilized using ethylene oxide		Caution
	Non-pyrogenic		Catalogue number		Unique Device Identifier
	Single sterile barrier system		Double sterile barrier system		Keep away from sunlight

Signs	Description	Signs	Description	Signs	Description
	Do not reuse		Do not use if the package is damaged		Keep dry
	Do not resterilize		Authorized representative in the European Community		Manufacturer
	CE marking of conformity, and Notified Body Code		Consult instructions for use	/	/

Manufacturer

Name: Haian Medigauze Co.,Ltd.

Address: Jiaoxie Industry park,226633Haian, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-513-88240958

Web: www.medigauze.com

E-mail: produce@medigauze.com

SRN: CN-MF-000009960

European Authorized Representative

Name: Shanghai International Holding Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537 Hamburg, Germany

E-mail: shholding@hotmail.com

SRN: DE-AR-000000001