



TIE-KTJKB-LKJ Assembly Instructions for Massage Bed


Parts List



Step 1: Install bed legs (7) onto the bed frame (8), (Note bottom screw direction) Fixate all 4 legs completely.



Step 2: Install tray (4) using screws (2), fixate shelf screws completely using lever (3). Note 3 horizontal pipes should be aligned with bottom rack.



Step 3: Stand the bed into upright direction, lift bed up and push down with force to secure bed legs on each side of the bed. Install backrest rod (1) into bed and frame sockets respectively.



Step 4: Fixate base tray onto horizontal rods (5), and place head pillow (6) onto the bed surface.



Recommended maximum load capacity for massage beds: 700lbs