

Electric Bed Installation Manual

Recommend to operate by 2 persons

Parts List



①



②



④



⑤



⑥



⑦



⑧



⑨



⑩



⑪



⑭

Prefixed Parts



③



⑫



⑬

Step 1: Place parts ⑨ and ⑩ as shown in the diagram, align the Velcro on the connecting parts, and press them together to secure them in place. Unscrew part ⑬ using part ⑭ screwdriver, turn the direction of the remote holder as shown below and fixate the part completely. As shown below:



Step 2: As shown in the diagram, use the provided iron rod and iron wire ⑫ to securely fasten the weighed iron base at the connection point between the motor and bed bracket by inserting the iron rod between the sockets and securing the iron wire at the end of the iron rod

As shown below:



Step 3: Prepare parts ① and ② as shown in the picture, totaling 8 sets of part ②. Then assemble the screws, spring pads, and flat pads in sequence before placing the bed frame on the mattress in the direction indicated in the picture. Ensure that the screw holes of the mattress and bed frame are completely fixed before securing them with set screws and tightening with part ①. Place part ⑪ in the position indicated in the diagram, press gently to secure it, and repeat the same operation for the other wheel.

As shown below:



Step 4: Prepare parts ③ and ④ and unscrew part ④ from the installed position. Insert part ③ into the screw, followed by the spring pad and flat pad in the same order. Place part ④ into the screw hole as shown in the picture, secure the screw, and tighten it using part ①.

As shown below:



Step 5: Align the junction of parts ⑤ and ⑦, pressing gently to ensure a firm bond. Next, insert the two brackets of part ⑦ into the two holes at the head of the bed, then secure them by pulling the buckle latch. As shown below:



Note: If you encounter any missing parts or other problems during installation, please feel free to contact us by email. Thank you for your support and purchase! Wish you a happy life!

Note: Recommended maximum load capacity for electric bed: 1000lbs