




D-22 Assembly Instructions

Parts List




Assembled by
Two Adults


Recommended maximum load capacity for massage beds: 660lbs

D-22 Assembly Instructions

Step 1: Insert the ⑩ ⑪ ⑫ accessory into the connecting piece of the ⑨ main frame as shown in the diagram, and tighten the ⑩ ⑪ square tube with the ⑬ short screw (Note: The middle of the 3 holes on each side of the ⑩ square tube is a through hole, and no screws need to be tightened).



Step 2: Connect the main frame of the other storage rack and tighten the screws.



Recommended maximum load capacity for massage beds: 660lbs

D-22 Assembly Instructions

Step 3: Use the ⑭ screwdriver to remove the hook of the ⑮ handle, rotate the direction of the hook, and install it as shown in the diagram. Use a large hexagonal wrench to remove the short screw at the end of the bed ⑯ and install the bed sheet stand ⑦. Install the backrest pneumatic lift switch onto the bed frame with ⑬ short screws.



Step 4: Install the ④ bed leg onto the bed frame with the ⑬ short screw (note the direction of the nut on the bed leg). Do not tighten the screw for the time being.




Recommended maximum load capacity for massage beds: 660lbs

D-22 Assembly Instructions

Step 5: Use ⑬ long screws to install the assembled shelf onto the bed legs, tighten the screws, and then tighten the screws used in step 4 to install the bed legs.



Step 6: Place the ⑧ prone hand board onto the ⑥ C-shaped pillow support and insert it into the holes at the head of the bed (the pull ropes on both sides are directly placed onto the aluminum rods of the C-shaped pillow support as shown in Figure 1, and the pull rope in the middle is installed onto the small round rods in the middle of the support as shown in Figures 2, 3, and 4).



1



2



3



4

Recommended maximum load capacity for massage beds: 660lbs

D-22 Assembly Instructions

Step 7: Place ⑤ C-shaped pillow on the C-shaped pillow support, ③ headrest on the bed surface, and ② storage board on the storage rack.



Recommended maximum load capacity for massage beds: 660lbs